

---

**readlif**

***Release 0.6.1***

**Apr 20, 2021**



---

## Contents:

---

<b>1 Indices and tables</b>	<b>5</b>
-----------------------------	----------

<b>Index</b>	<b>7</b>
--------------	----------



---

```
class readlif.reader.LifFile(filename)
```

Given a path or buffer to a lif file, returns objects containing the image and data.

This is based on the java openmicroscopy bioformats lif reading code that is here: <https://github.com/openmicroscopy/bioformats/blob/master/components/formats-gpl/src/loci/formats/in/LIFReader.java> # noqa

### Variables

- **xml\_header** (*string*) – The LIF xml header with tons of data
- **xml\_root** (*ElementTree*) – ElementTree XML representation
- **offsets** (*list*) – Byte positions of the files
- **num\_images** (*int*) – Number of images
- **image\_list** (*dict*) – Has the keys: path, folder\_name, folder\_uuid, name, image\_id, frames

### Example

```
>>> from readlif.reader import LifFile
>>> new = LifFile('./path/to/file.lif')
```

```
>>> for image in new.get_iter_image():
>>>     for frame in image.get_iter_t():
>>>         frame.image_info['name']
>>>         # do stuff
```

```
>>> # For non-xy imaging experiments
>>> img_0 = new.get_image(0)
>>> for i in range(0, img_0.dims_n[4]):
>>>     plane = img_0.get_plane(requested_dims = {4: i}))
```

#### **get\_image** (*img\_n=0*)

Specify the image number, and this returns a LifImage object of that image.

**Parameters** **img\_n** (*int*) – Image number to retrieve

**Returns** LifImage object with specified image

#### **get\_iter\_image** (*img\_n=0*)

Returns an iterator of LifImage objects in the lif file.

**Parameters** **img\_n** (*int*) – Image to start iteration at

**Returns** Iterator of LifImage objects.

```
class readlif.reader.LifImage(image_info, offsets, filename)
```

This should not be called directly. This should be generated while calling get\_image or get\_iter\_image from a LifFile object.

### Variables

- **path** (*str*) – path / name of the image
- **dims** (*tuple*) – (x, y, z, t, m)
- **display\_dims** (*tuple*) – The first two dimensions of the lif file. This is used to decide what dimensions are returned in a 2D plane.

- **dims\_n** (*dict*) – {0: length, 1: length, 2: length, n: length}

For atypical imaging experiments, i.e. those not simple photos of XY frames, this attribute will be more useful than *dims*. This attribute will hold a dictionary with the length of each dimension, in the order it is referenced in the .lif file.

Currently, only some of the 10 possible dimensions are used / known:

- 1: x-axis
- 2: y-axis
- 3: z-axis
- 4: time
- 5: detection wavelength
- 6-8: Unknown
- 9: illumination wavelength
- 10: mosaic tile

- **name** (*str*) – image name

- **offsets** (*list*) – Byte position offsets for each image.

- **filename** (*str, bytes, os.PathLike, io.IOBase*) – The name of the LIF file being read

- **channels** (*int*) – Number of channels in the image

- **nz** (*int*) – number of ‘z’ frames

Note, it is recommended to use *dims.z* instead. However, this will be kept for compatibility.

- **nt** (*int*) – number of ‘t’ frames

Note, it is recommended to use *dims.t* instead. However, this will be kept for compatibility.

- **scale** (*tuple*) – (scale\_x, scale\_y, scale\_z, scale\_t).

Conversion factor: px/µm for x, y and z; sec/image for t.

- **scale\_n** (*dict*) – {0: length, 1: length, 2: length...}.

Conversion factor: px/µm for x, y and z; sec/image for t. Related to *dims\_n* above.

- **bit\_depth** (*tuple*) – A tuple of ints that indicates the bit depth of each channel in the image.

- **mosaic\_position** (*list*) – If the image is a mosaic (tiled), this contains a list of tuples with four values: (*FieldX, FieldY, PosX, PosY*). The length of this list is equal to the number of tiles.

- **info** (*dict*) – Direct access to data dict from LifFile, this is most useful for debugging. These are values pulled from the Leica XML.

### **get\_frame** (*z=0, t=0, c=0, m=0*)

Gets the specified frame (z, t, c, m) from image.

#### **Parameters**

- **z** (*int*) – z position
- **t** (*int*) – time point
- **c** (*int*) – channel

- **m** (*int*) – mosaic image

**Returns** Pillow Image object

**get\_iter\_c** (*z=0, t=0, m=0*)

Returns an iterator over the channels at time t and position z.

**Parameters**

- **z** (*int*) – z position
- **t** (*int*) – time point
- **m** (*int*) – mosaic image

**Returns** Iterator of Pillow Image objects

**get\_iter\_m** (*z=0, t=0, c=0*)

Returns an iterator over the z series of time t and channel c.

**Parameters**

- **t** (*int*) – time point
- **c** (*int*) – channel
- **z** (*int*) – z position

**Returns** Iterator of Pillow Image objects

**get\_iter\_t** (*z=0, c=0, m=0*)

Returns an iterator over time t at position z and channel c.

**Parameters**

- **z** (*int*) – z position
- **c** (*int*) – channel
- **m** (*int*) – mosaic image

**Returns** Iterator of Pillow Image objects

**get\_iter\_z** (*t=0, c=0, m=0*)

Returns an iterator over the z series of time t and channel c.

**Parameters**

- **t** (*int*) – time point
- **c** (*int*) – channel
- **m** (*int*) – mosaic image

**Returns** Iterator of Pillow Image objects

**get\_plane** (*display\_dims=None, c=0, requested\_dims=None*)

Gets the specified frame from image.

**Parameters**

- **display\_dims** (*tuple*) – Two value tuple (1, 2) specifying the two dimension plane to return. This will default to the first two dimensions in the LifFile, specified by LifFile.display\_dims
- **c** (*int*) – channel
- **requested\_dims** (*dict*) – Dictionary indicating the item to be returned, as described by a numerical dictionary, ex: {3: 0, 4: 1}

**Returns** Pillow Image object

`readlif.utilities.get_xml (filename)`

Given a lif file, returns two values (xml\_root, xml\_header) where xml\_root is an ElementTree root, and xml\_header is the text.

This is useful for debugging.

Some private functions are used from readlif.reader.

**Parameters** `filename` (*string*) – what file to open?

# CHAPTER 1

---

## Indices and tables

---

- genindex
- modindex
- search



---

## Index

---

### G

`get_frame()` (*readlif.reader.LifImage method*), 2  
`get_image()` (*readlif.reader.LifFile method*), 1  
`get_iter_c()` (*readlif.reader.LifImage method*), 3  
`get_iter_image()` (*readlif.reader.LifFile method*), 1  
`get_iter_m()` (*readlif.reader.LifImage method*), 3  
`get_iter_t()` (*readlif.reader.LifImage method*), 3  
`get_iter_z()` (*readlif.reader.LifImage method*), 3  
`get_plane()` (*readlif.reader.LifImage method*), 3  
`get_xml()` (*in module readlif.utilities*), 4

### L

`LifFile` (*class in readlif.reader*), 1  
`LifImage` (*class in readlif.reader*), 1

### R

`readlif.reader` (*module*), 1  
`readlif.utilities` (*module*), 4